

Quality right down to the smallest detail

Liebherr Mediline laboratory refrigerators the ideal solution for a small footprint or for under-worktop integration. The range comprises two freestanding and under-worktop refrigerators with a glass door or a solid door. The precision Comfort electronic

controller enables temperatures to be set accurately. The dynamic cooling system ensures high temperature consistency. Integrated alarm systems guarantee safe storage.



The Comfort electronic controller.

The precision Comfort electronic controller with digital temperature display enables temperatures to be set accurately. Symbols indicate the operating status of the appliance. The electronic controller is flush-mounted and equipped with a membrane keyboard for perfect hygiene.




Integrated alarm systems.

Visual and audible alarm systems alert the user to undesired temperature deviations or when the door is open. All alarm parameters can be individually set. For example, the door-open alarm delay can be adjusted between 1 and 5 minutes. Additional safety is provided by a visual power failure alarm and sensor defect alarm.




External temperature and alarm documentation.

The laboratory appliances are equipped with a volt-free contact for alarm forwarding to an external remote warning system, as well as with a serial interface RS 485 which enables central documentation of temperature profiles and alarm events. LTM documentation software is optionally available.



Maximum temperature consistency.

The dynamic cooling system in conjunction with the precision Comfort electronic controller ensures maximum temperature stability. All laboratory appliances with a Comfort electronic controller are designed according to the ISO 60068-3 standard with regard to optimum temperature consistency.




Accurate 1-point calibration.

The laboratory appliances with a Comfort electronic controller are equipped with 1-point calibration for accurate temperature control. This permits compensation between the set temperature and the actual internal temperature. Positive or negative changes to the correction value can be made in steps of 0.1 K.



Integrated data memory.

The Comfort electronic controller is equipped with an integrated data memory which logs the maximum and minimum internal temperatures as well as the last three temperature alarms and power failures, together with the time, date and duration. These data can be read out on the display.



Safety thermostat.

The laboratory appliances are equipped with an additional safety thermostat to prevent temperatures below +2°C and to protect sensitive products in the event of a fault.



External temperature sensor.

The laboratory appliances with a Comfort electronic controller have an access port (diameter 10 mm) on the rear which allows the temperature sensor to be positioned in the interior.