

Refrigeration unit

PAS 10001



Independent refrigeration unit

PAS10001

This independent refrigeration unit has been designed and developed to meet standards on human milk cooling.

This equipment should be connected to a HSC pasteuriser PAS10000 model. It runs automatically when the pasteurisation holding time is over. It operates a full cooling cycle down to 4 °C to meet actual regulation on human milk treatment and preserve the inactivation of pathogens and avoid re-growth of bacterial flora inside pasteurised human milk.



The pasteuriser PAS10000 model associated with the refrigeration unit becomes the PAS10001 model.

- Integrated water reserve
- Cooling duration at 4 °C: less than one hour
- Full automatic operation
- Follow up of the cooling cycle on the touch screen interface
- Adjustable cooling temperature

